BLASTX™

Description

Next Science BlastX is an antimicrobial wound gel powered by a non-toxic technology that inhibits the growth of microorganisms within the gel application.

BlastX contains Benzalkonium Chloride 0.13%, Polyethylene Glycol 400, Polyethylene Glycol 3350, Sodium Citrate, Citric Acid, and Water.

Indications for Use

BlastX is indicated for the management of wounds such as stage I-IV pressure ulcers, partial- and full-thickness wounds, diabetic foot and leg ulcers, post-surgical wounds, first and second degree burns, and grafted and donor sites.

Contraindications

BlastX should not be used if there is a history of allergy to any of the ingredients.

Warnings and Precautions

For external topical use only. Do not ingest. Must not be used in or near the eyes. Do not cover with alginate dressing.

Instructions of Use

BlastX should be applied directly to the wound and then cover with an appropriate dressing as needed for up to 14 uses within a 30 day period. After opening, product may be used until labeled expiration date.


www.nextscience.com/wgpatents

Powered by X^{BIO} Technology

Manufactured for:
Next Science, LLC
10550 Deerwood Park Blvd
Suite 300
Jacksonville, FL 32256 USA

www.nextscience.com

15°C  25°C
Rx Only

 Consult instructions for use * Reg #1641	 Manufacturer * Reg #3082
 Use-by-date * Reg #2607	 Temperature storage limit * Reg #0632
 Batch code * Reg #2492	 Caution * Reg #0434A
 Catalogue Number * Reg #2493	 Net Weight
Rx Only	CAUTION: Federal (USA) law restricts this product to sale and use by, or on the order of, a physician.

* ISO 15223-1 - Medical Devices - symbols to be used with medical device labels, labeling and information to be supplied

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